



Mouse Anti-763.74 monoclonal antibody, clone MK2-23 (CABT-ZB1185)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Monoclonal antibody 763.74.
Isotype	IgG1, κ
Source/Host	Mouse
Clone	MK2-23
Purification	Protein A purified
Conjugate	Unconjugated
Format	Liquid
Concentration	Lot specific
Size	200 μ g
Buffer	PBS
Preservative	0.02% Proclin 300
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Ship	Wet ice