



Human Anti-Trastuzumab monoclonal antibody, clone scFv69 (CABT-ZB1186)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Trastuzumab F(ab') ₂ fragments.
Isotype	IgG1, λ
Source/Host	Human
Clone	scFv69
Purification	Protein A purified
Conjugate	Unconjugated
Applications	Suitable for Blocking. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	200 µg
Buffer	PBS
Preservative	0.02% Proclin 300
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Ship	Wet ice