



# Mouse Anti-Gemtuzumab monoclonal antibody, clone 975155 (CABT-ZB1190)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	No cross-reactivity with Cetuximab, Ibritumomab, Rituximab, Tositumomab, Brentuximab, Trastuzumab, Solanezumab, Teprotumumab, and Aducanumab was observed.
<b>Immunogen</b>	Full length Gemtuzumab
<b>Isotype</b>	IgG2a
<b>Source/Host</b>	Mouse
<b>Clone</b>	975155
<b>Purification</b>	Protein A or G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Suitable for WB, ELISA, FC. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Reconstitution</b>	Reconstitute at 0.8 mg/mL in sterile PBS.
<b>Format</b>	Lyophilized
<b>Size</b>	100 µg
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	12 months from date of receipt, -20 to -70 °C as supplied.

1 month, 2 to 8 °C under sterile conditions after reconstitution.  
6 months, -20 to -70 °C under sterile conditions after reconstitution.

---

**Ship**

Wet ice

---