



Mouse Anti-Gemtuzumab monoclonal antibody, clone 975155 [Alexa Fluor 488] (CABT-ZB1191)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	No cross-reactivity with Cetuximab, Ibritumomab, Rituximab, Tositumomab, Brentuximab, Trastuzumab, Solanezumab, Teprotumumab, and Aducanumab was observed.
Immunogen	Full length Gemtuzumab
Isotype	IgG2a
Source/Host	Mouse
Clone	975155
Purification	Protein A or G purified
Conjugate	Alexa Fluor 488
Applications	Recommended Working Concentration: WB: 0.5 µg/mL FC: 0.25 µg/10 ⁶ cells Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Reconstitution	Reconstitute at 0.9 mg/mL in sterile PBS.
Format	Lyophilized
Size	100 µg
Buffer	PBS

Preservative	None
Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied.
Ship	Wet ice