



Mouse Anti-Ibritumomab monoclonal antibody, clone 974622 (CABT-ZB1192)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	100% cross-reactivity with Tositumomab and no cross-reactivity with Cetuximab, Rituximab, Trastuzumab, Brentuximab, Gemtuzumab, Solanezumab, Teprotumumab, and Aducanumab was observed.
Immunogen	Full length Ibritumomab
Isotype	IgG2a
Source/Host	Mouse
Clone	974622
Purification	Protein A or G purified
Conjugate	Unconjugated
Applications	Suitable for WB, ELISA, FC. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Reconstitution	Reconstitute at 0.10 mg/mL in sterile PBS.
Format	Lyophilized
Size	100 µg
Buffer	PBS
Preservative	None

Storage	12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution.
----------------	---

Ship	Wet ice
-------------	---------
