



Rabbit Anti-Solanezumab monoclonal antibody, clone 2372A (CABT-ZB1194)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length Solanezumab
Isotype	IgG
Source/Host	Rabbit
Clone	2372A
Purification	Protein A or G purified
Conjugate	Unconjugated
Applications	Suitable for WB, ELISA, FC. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Reconstitution	Reconstitute at 0.12 mg/mL in sterile PBS.
Format	Lyophilized
Size	100 µg
Buffer	PBS
Preservative	None
Storage	12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution.

Ship

Wet ice
