



Rabbit Anti-Solanezumab monoclonal antibody, clone 2372A [Alexa Fluor 488] (CABT-ZB1195)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length Solanezumab
Isotype	IgG
Source/Host	Rabbit
Clone	2372A
Purification	Protein A or G purified
Conjugate	Alexa Fluor 488
Applications	Recommended Working Concentration: WB: 1 µg/mL FC: 0.25 µg/10 ⁶ cells Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Reconstitution	Reconstitute at 0.13 mg/mL in sterile PBS.
Format	Lyophilized
Size	100 µg
Buffer	PBS
Preservative	None

Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied.
Ship	Wet ice