



Rabbit Anti-Teplizumab monoclonal antibody, clone 2505B [Alexa Fluor 488] (CABT-ZB1197)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length Teplizumab
Isotype	IgG
Source/Host	Rabbit
Clone	2505B
Purification	Protein A or G purified
Conjugate	Alexa Fluor 488
Applications	Recommended Working Concentration: WB: 2 µg/mL FC: 0.25 µg/10 ⁶ cells Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Reconstitution	Reconstitute at 0.15 mg/mL in sterile PBS.
Format	Lyophilized
Size	100 µg
Buffer	PBS
Preservative	None
Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied.

Ship

Wet ice
