



# Mouse Anti-Teprotumumab monoclonal antibody, clone 975547 (CABT-ZB1198)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Full length Teprotumumab
<b>Isotype</b>	IgG2a
<b>Source/Host</b>	Mouse
<b>Clone</b>	975547
<b>Purification</b>	Protein A or G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Suitable for ELISA, FC. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Reconstitution</b>	Reconstitute at 0.16 mg/mL in sterile PBS.
<b>Format</b>	Lyophilized
<b>Size</b>	100 µg
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution.

Ship

Wet ice

---