



Chicken Anti-Rituximab polyclonal antibody (CABT-ZB1202)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgY
Source/Host	Chicken
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	PBS, with 0.02% sodium azide, pH7
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for up to 1 month. For longer storage, aliquot and store at -20°C.
Ship	Wet ice