



Chicken Anti-Tocilizumab polyclonal antibody (CABT-ZB1203)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgY
Source/Host	Chicken
Conjugate	Unconjugated
Applications	Suitable for ELISA, Blocking Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	PBS, with 0.02% sodium azide, pH7
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for up to 1 month. For longer storage, aliquot and store at -20°C.
Ship	Wet ice