



Rabbit Anti-Human ZDHC8 polyclonal antibody (CABT-ZB121)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects endogenous levels of ZDHC8 at Human, Mouse
Target	ZDHC8
Immunogen	Synthesized peptide derived from human ZDHC8
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB Recommended Dilution: WB: 1:500-2000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS with 50% Glycerol

Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	ZDHHC8; zinc finger, DHHC-type containing 8; probable palmitoyltransferase ZDHHC8; KIAA1292; ZNF378; zinc finger protein 378; zinc finger, DHHC domain like containing 1; membrane-associated DHHC8 zinc finger protein; DHHC8; ZDHHCL1
-----------------	---

GENE INFORMATION

Synonyms	ZDHHC8; zinc finger, DHHC-type containing 8; probable palmitoyltransferase ZDHHC8; KIAA1292; ZNF378; zinc finger protein 378; zinc finger, DHHC domain like containing 1; membrane-associated DHHC8 zinc finger protein; DHHC8; ZDHHCL1
-----------------	---

Entrez Gene ID	29801
-----------------------	-----------------------

UniProt ID	Q9ULC8
-------------------	------------------------