



Rabbit Anti-Human ZDHC20 (aa 72-121) polyclonal antibody (CABT-ZB122)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects endogenous levels of ZDHC20 at Human, Mouse, Rat
Target	ZDHC20
Immunogen	Synthesized peptide derived from human ZDHC20 (aa 72-121).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA Recommended Dilution: WB: 1:500-2000 ELISA: 1:10000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS with 50% Glycerol

Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	ZDHHC20; zinc finger, DHHC-type containing 20; probable palmitoyltransferase ZDHHC20; FLJ25952; DHHC-20; 4933421L13Rik; DHHC-containing protein 20; zinc finger, DHHC domain containing 20; zinc finger DHHC domain-containing protein 20; MGC126005
-----------------	--

GENE INFORMATION

Synonyms	ZDHHC20; zinc finger, DHHC-type containing 20; probable palmitoyltransferase ZDHHC20; FLJ25952; DHHC-20; 4933421L13Rik; DHHC-containing protein 20; zinc finger, DHHC domain containing 20; zinc finger DHHC domain-containing protein 20; MGC126005
Entrez Gene ID	253832
UniProt ID	Q5W0Z9