



# Rabbit Anti-Human ZDHC18 (aa 20-100) polyclonal antibody (CABT-ZB128)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	This antibody detects endogenous levels of ZDHC18 at Human, Mouse, Rat
<b>Target</b>	ZDHC18
<b>Immunogen</b>	Synthesized peptide derived from human ZDHC18 (aa 20-100).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA Recommended Dilution: WB: 1:500-2000 ELISA: 1:5000-20000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	100 µg
<b>Buffer</b>	PBS with 50% Glycerol

<b>Preservative</b>	0.02% sodium azide
<b>Storage</b>	Long time storage is recommended at -20°C.
<b>Ship</b>	Wet ice

## BACKGROUND

<b>Keywords</b>	ZDHHC18; zinc finger, DHHC-type containing 18; palmitoyltransferase ZDHHC18; DKFZp667O2416; zinc finger, DHHC domain containing 18; zinc finger DHHC domain-containing protein 18; DHHC18; DHHC-18
-----------------	--

## GENE INFORMATION

<b>Synonyms</b>	ZDHHC18; zinc finger, DHHC-type containing 18; palmitoyltransferase ZDHHC18; DKFZp667O2416; zinc finger, DHHC domain containing 18; zinc finger DHHC domain-containing protein 18; DHHC18; DHHC-18
-----------------	--

<b>Entrez Gene ID</b>	<a href="#">84243</a>
-----------------------	-----------------------

<b>UniProt ID</b>	<a href="#">Q9NUE0</a>
-------------------	------------------------