



Rabbit Anti-Human LYPLA2 (aa 30-110) polyclonal antibody (CABT-ZB130)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects endogenous levels of LYPLA2 at Human, Mouse, Rat
Target	LYPLA2
Immunogen	Synthesized peptide derived from human LYPLA2 (aa 30-110).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA Recommended Dilution: WB: 1:500-2000 ELISA: 1:5000-20000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS with 50% Glycerol

Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	LYPLA2; lysophospholipase II; acyl-protein thioesterase 2; APT 2; LPL-II; lysoPLA II; APT-2; DJ886K2.4
-----------------	--

GENE INFORMATION

Synonyms	LYPLA2; lysophospholipase II; acyl-protein thioesterase 2; APT 2; LPL-II; lysoPLA II; APT-2; DJ886K2.4
-----------------	--

Entrez Gene ID	11313
-----------------------	-----------------------

UniProt ID	O95372
-------------------	------------------------