



# Rabbit Anti-Human ABHD8 (aa 291-340) polyclonal antibody (CABT-ZB133)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	This antibody detects endogenous levels of ABHD8 at Human, Mouse
<b>Target</b>	ABHD8
<b>Immunogen</b>	Synthesized peptide derived from human ABHD8 (aa 291-340).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC-P, ELISA Recommended Dilution: WB: 1:500-2000 IHC: 1:100 - 1:300 ELISA: 1:40000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	100 µg

<b>Buffer</b>	PBS with 50% Glycerol
<b>Preservative</b>	0.02% sodium azide
<b>Storage</b>	Long time storage is recommended at -20°C.
<b>Ship</b>	Wet ice

## BACKGROUND

**Keywords** ABHD8; Abhydrolase Domain Containing 8; Alpha/Beta Hydrolase Domain-Containing Protein 8; Abhydrolase Domain-Containing Protein 8; Protein ABHD8; FLJ11743; MGC14280; MGC2512

## GENE INFORMATION

**Synonyms** ABHD8; Abhydrolase Domain Containing 8; Alpha/Beta Hydrolase Domain-Containing Protein 8; Abhydrolase Domain-Containing Protein 8; Protein ABHD8; FLJ11743; MGC14280; MGC2512

**Entrez Gene ID** [79575](#)

**UniProt ID** [Q96113](#)