



Rabbit Anti-Human ABHD2 (aa 251-300) polyclonal antibody (CABT-ZB136)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects endogenous levels of ABHD2 at Human, Mouse
Target	ABHD2
Immunogen	Synthesized peptide derived from human ABHD2 (aa 251-300).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IF, ELISA Recommended Dilution: WB: 1:500-2000 IF: 1:200-1:1000 ELISA: 1:40000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg

Buffer	PBS with 50% Glycerol
Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	ABHD2; abhydrolase domain containing 2; abhydrolase domain-containing protein 2; LABH2; protein PHPS1-2; lung alpha/beta hydrolase 2; alpha/beta hydrolase domain containing protein 2; HS1-2; PHPS1-2; MGC26249; MGC111112
-----------------	---

GENE INFORMATION

Synonyms	ABHD2; abhydrolase domain containing 2; abhydrolase domain-containing protein 2; LABH2; protein PHPS1-2; lung alpha/beta hydrolase 2; alpha/beta hydrolase domain containing protein 2; HS1-2; PHPS1-2; MGC26249; MGC111112
-----------------	---

Entrez Gene ID	11057
-----------------------	-----------------------

UniProt ID	P08910
-------------------	------------------------