



Rabbit Anti-Human ABHD12B (aa 281-330) polyclonal antibody (CABT-ZB138)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects endogenous levels of ABHD12B at Human
Target	ABHD12B
Immunogen	Synthesized peptide derived from human ABHD12B (aa 281-330).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purification
Conjugate	Unconjugated
Applications	IHC-P, ELISA Recommended Dilution: IHC: 1:100 - 1:300 ELISA: 1:40000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS with 50% Glycerol

Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	ABHD12B; abhydrolase domain containing 12B; BEM46L3; C14orf29; c14_5314; abhydrolase domain-containing protein 12B
-----------------	--

GENE INFORMATION

Synonyms	ABHD12B; abhydrolase domain containing 12B; BEM46L3; C14orf29; c14_5314; abhydrolase domain-containing protein 12B
-----------------	--

Entrez Gene ID	145447
-----------------------	------------------------

UniProt ID	Q7Z5M8
-------------------	------------------------