



Rabbit Anti-Human ABHD3 (aa 240-320) polyclonal antibody (CABT-ZB141)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects endogenous levels of ABHD3 at Human, Mouse
Target	ABHD3
Immunogen	Synthesized peptide derived from human ABHD3 (aa 240-320).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA Recommended Dilution: WB: 1:500-2000 ELISA: 1:5000-20000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS with 50% Glycerol

Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	ABHD3; Abhydrolase Domain Containing 3, Phospholipase; LABH3; Abhydrolase Domain-Containing Protein; Phospholipase ABHD3; Alpha/Beta Hydrolase Domain Containing Protein 3; Abhydrolase Domain Containing 3; Lung Alpha/Beta Hydrolase 3
-----------------	--

GENE INFORMATION

Synonyms	ABHD3; Abhydrolase Domain Containing 3, Phospholipase; LABH3; Abhydrolase Domain-Containing Protein; Phospholipase ABHD3; Alpha/Beta Hydrolase Domain Containing Protein 3; Abhydrolase Domain Containing 3; Lung Alpha/Beta Hydrolase 3
Entrez Gene ID	171586
UniProt ID	Q8WU67