



Rabbit Anti-Human PPT2 polyclonal antibody (CABT-ZB145)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody detects endogenous levels of PPT2 at Human, Mouse, Rat
Target	PPT2
Immunogen	Synthesized peptide derived from human PPT2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB Recommended Dilution: WB: 1:500-2000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS with 50% Glycerol

Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	PPT2; palmitoyl-protein thioesterase 2; lysosomal thioesterase PPT2; VE-statin2; S-thioesterase G14; palmitoyl-protein hydrolase 2; G14; PPT-2; C6orf8; DKFZp564P1516
-----------------	---

GENE INFORMATION

Synonyms	PPT2; palmitoyl-protein thioesterase 2; lysosomal thioesterase PPT2; VE-statin2; S-thioesterase G14; palmitoyl-protein hydrolase 2; G14; PPT-2; C6orf8; DKFZp564P1516
Entrez Gene ID	9374
UniProt ID	Q9UMR5