



Human Anti-Human TK1 monoclonal antibody, Clone A11 (CABT-ZB146)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	TK1
Isotype	IgG
Source/Host	Human
Species Reactivity	Human
Clone	A11
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA Capture We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB146 - CABT-ZB147; CABT-ZB146 - CABT-ZB148; CABT-ZB146 - CABT-ZB149
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords TK1; thymidine kinase 1, soluble; TK2; thymidine kinase, cytosolic; thymidine kinase-1;
thymidine kinase 1 soluble isoform

GENE INFORMATION

Synonyms TK1; thymidine kinase 1, soluble; TK2; thymidine kinase, cytosolic; thymidine kinase-1;
thymidine kinase 1 soluble isoform

Entrez Gene ID [7083](#)

UniProt ID [P04183](#)
