



Human Anti-Human TK1 monoclonal antibody, Clone B22 (CABT-ZB147)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	TK1
Isotype	IgG
Source/Host	Human
Species Reactivity	Human
Clone	B22
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA Detection We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB146 - CABT-ZB147
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Keywords	TK1; thymidine kinase 1, soluble; TK2; thymidine kinase, cytosolic; thymidine kinase-1; thymidine kinase 1 soluble isoform
-----------------	--

GENE INFORMATION

Synonyms	TK1; thymidine kinase 1, soluble; TK2; thymidine kinase, cytosolic; thymidine kinase-1; thymidine kinase 1 soluble isoform
-----------------	--

Entrez Gene ID	7083
-----------------------	----------------------

UniProt ID	P04183
-------------------	------------------------
