



Sheep Anti-17 α -Boldenone polyclonal antibody (CABT-ZB151)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Boldenone
Immunogen	17 α -Boldenone-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	20 mM Phosphate, 150 mM Sodium Chloride.
Preservative	0.09% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Boldenone is an anabolic steroid developed for veterinary use, mostly for treatment of horses. It is not indicated for use in humans in the US and is only available through veterinary clinics.
Keywords	Boldenone; Dehydrotestosterone; 846-48-0; 1-Dehydrotestosterone; 17beta-Boldenone

GENE INFORMATION

Synonyms	Boldenone; Dehydrotestosterone; 846-48-0; 1-Dehydrotestosterone; 17beta-Boldenone
-----------------	---