



# Sheep Anti-3-Acetyl-deoxynivalenol polyclonal antibody (CABT-ZB154)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	Deoxynivalenol
<b>Immunogen</b>	3-Acetyl-DON-KLH
<b>Isotype</b>	IgG
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	N/A
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg
<b>Buffer</b>	20 mM Phosphate, 150 mM Sodium Chloride.
<b>Preservative</b>	0.09% sodium azide
<b>Storage</b>	Long time storage is recommended at -20°C.
<b>Ship</b>	Wet ice

## BACKGROUND

**Introduction**

3-Acetyl Deoxynivalenol is a natural product found in *Fusarium*, *Fusarium culmorum*, and *Fusarium graminearum* with data available.

---

**Keywords**

3-Acetyl Deoxynivalenol

---

## GENE INFORMATION

**Synonyms**

3-Acetyl Deoxynivalenol

---