



# Sheep Anti-Dihydrotestosterone polyclonal antibody (CABT-ZB156)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	DHT
Immunogen	Dihydrotestosterone-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	20 mM Phosphate, 150 mM Sodium Chloride.
Preservative	0.09% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

## BACKGROUND

**Introduction**

A potent androgenic metabolite of testosterone. Dihydrotestosterone (DHT) is generated by a 5-alpha reduction of testosterone. Unlike testosterone, DHT cannot be aromatized to estradiol therefore DHT is considered a pure androgenic steroid.

**Keywords**

DHT; Dihydrotestosterone

## GENE INFORMATION

**Synonyms**

DHT; Dihydrotestosterone