



Sheep Anti-Epitestosterone polyclonal antibody (CABT-ZB157)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Epitestosterone
Immunogen	Epi-testosterone-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	20 mM Phosphate, 150 mM Sodium Chloride.
Preservative	0.09% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Epitestosterone is an androstanoid that is the C-17 epimer of testosterone. It has a role as an androgen antagonist and a human metabolite. It is an androstanoid, a 17alpha-hydroxy steroid and a 3-oxo-Delta(4) steroid.

Keywords

Epitestosterone; 481-30-1; Isotestosterone; cis-Testosterone; 17-Epitestosterone

GENE INFORMATION

Synonyms

Epitestosterone; 481-30-1; Isotestosterone; cis-Testosterone; 17-Epitestosterone