



Sheep Anti-Methylboldenone/Methyltestosterone polyclonal antibody (CABT-ZB158)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Methyltestosterone
Immunogen	Medroxyprogesterone Acetate-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	20 mM Phosphate, 150 mM Sodium Chloride.
Preservative	0.09% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Methyltestosterone is a 17beta-hydroxy steroid that is testosterone bearing a methyl group at the 17alpha position. It has a role as an antineoplastic agent, an anabolic agent and an androgen. It is a 3-oxo-Delta(4) steroid, a 17beta-hydroxy steroid and an enone. It derives from a testosterone.
Keywords	methyltestosterone; 17-Methyltestosterone; Testred; Metandren; Android

GENE INFORMATION

Synonyms	methyltestosterone; 17-Methyltestosterone; Testred; Metandren; Android
----------	--