



# Sheep Anti-Ochratoxin A polyclonal antibody (CABT-ZB159)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	Ochratoxin A
<b>Immunogen</b>	Ochratoxin A-KLH
<b>Isotype</b>	IgG
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	N/A
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg
<b>Buffer</b>	20 mM Phosphate, 150 mM Sodium Chloride.
<b>Preservative</b>	0.09% sodium azide
<b>Storage</b>	Long time storage is recommended at -20°C.
<b>Ship</b>	Wet ice

## BACKGROUND

**Introduction** Ochratoxin A—a toxin produced by different *Aspergillus* and *Penicillium* species — is one of the most-abundant food-contaminating mycotoxins.

---

**Keywords** OCHRATOXIN A; 303-47-9; Phenylalanine - ochratoxin A; Antibiotic 9663; 3R14S-Ochratoxin A

---

## GENE INFORMATION

**Synonyms** OCHRATOXIN A; 303-47-9; Phenylalanine - ochratoxin A; Antibiotic 9663; 3R14S-Ochratoxin A

---