



Sheep Anti-Tilidine polyclonal antibody (CABT-ZB162)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Tilidine
Immunogen	Tilidine-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	20 mM Phosphate, 150 mM Sodium Chloride.
Preservative	0.09% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Tilidine, or tilidate (brand names: Tilidin, Valoron and Valtran) is a synthetic opioid painkiller, used mainly in Belgium, Bulgaria, Germany, Luxembourg, South Africa and Switzerland for the treatment of moderate to severe pain, both acute and chronic.

Keywords

Dextilidine; TILIDINE; Dextilidina; Dextilidinum; 32447-90-8

GENE INFORMATION

Synonyms

Dextilidine; TILIDINE; Dextilidina; Dextilidinum; 32447-90-8
