



Sheep Anti-Secobarbital monoclonal antibody, clone 21309 (CABT-ZB165)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Secobarbital
Immunogen	Secobarbital-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Clone	21309
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	20 mM Phosphate, 150 mM Sodium Chloride.
Preservative	0.09% sodium azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Secobarbital (marketed by Eli Lilly and Company under the brand names Seconal and Tuinal) is a barbiturate derivative drug with anaesthetic, anticonvulsant, sedative and hypnotic properties. It is commonly known as quinalbarbitone in the United Kingdom.
Keywords	secobarbital; Quinalbarbitone; Secobarbitone; Seconal; Meballymal

GENE INFORMATION

Synonyms	secobarbital; Quinalbarbitone; Secobarbitone; Seconal; Meballymal
----------	---