



Goat Anti-Streptococcus Pyogenes polyclonal antibody (CABT-ZB279)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Streptococcus Pyogenes
Immunogen	Streptococcus Pyogenes
Isotype	IgG
Source/Host	Goat
Species Reactivity	Streptococcus
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IF
Reconstitution	Reconstitution by adding 1.0 mL of reagent quality water to the product vial.
Format	Purified, Lyophilized
Size	1 mg
Buffer	PBS
Preservative	None
Storage	<p>Lyophilized material should be kept at 4°C until rehydrated and will remain stable for a minimum of 1 year.</p> <p>Store the reconstitution liquid solution at 4°C for up to a week; thereafter, aliquots should be prepared and stored frozen. Stable for a minimum of 1 year at -20°C.</p>

Ship

Wet ice

BACKGROUND

Introduction

Streptococcus is a genus of spherical Gram-positive bacteria belonging to the phylum Firmicutes and the lactic acid bacteria group.

Keywords

Streptococcal groups; Streptococcaceae; Streptococcus; GAS

GENE INFORMATION

Synonyms

Streptococcal groups; Streptococcaceae; Streptococcus; GAS
