



Mouse Anti-Human CD130 monoclonal antibody, clone C-U3 (CABT-ZB307)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes the Gp130, common subunit for IL-6, IL-11, OSM, LIF, CNTF, CT-1 receptors, a 130-140 kDa protein
Target	IL6ST
Immunogen	Natural soluble gp130
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-U3
Conjugate	Unconjugated
Applications	ELISA(Det) Recommended Working Concentration: ELISA: 0.05-0.5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB305 - CABT-ZB307
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS

Preservative	None
Storage	Stable at +2-8°C for 12 months
Ship	Wet ice

GENE INFORMATION

Synonyms	IL6ST; interleukin 6 signal transducer; CD130; GP130; CDW130; IL-6RB; interleukin-6 receptor subunit beta; CD130 antigen; IL-6R subunit beta; membrane glycoprotein 130
Entrez Gene ID	16195
UniProt ID	Q00560