



Mouse Anti-Human CD141 monoclonal antibody, clone C-B46 [FITC] (CABT-ZB311)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes thrombomodulin (TM), a 105 kDa protein
Target	THBD
Immunogen	TF1 cell line
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-B46
Conjugate	FITC
Applications	ELISA(Cap), FC Recommended Working Concentration: ELISA: 0.2-1 µg/ml FC: Use 10 µl to label 10⁶ cells or 100 µl of whole blood We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB309 - CABT-ZB313 CABT-ZB311 - CABT-ZB313
Format	Purified, Liquid
Concentration	Lot specific
Size	100 tests (1 mL)
Buffer	PBS

Preservative	0.1% Sodium Azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Ship	Wet ice

GENE INFORMATION

Synonyms	THBD; thrombomodulin; TM; THRM; AHUS6; BDCA3; CD141; THPH12; fetomodulin; CD141 antigen
Entrez Gene ID	7056
UniProt ID	P07204