



Mouse Anti-Human CD141 monoclonal antibody, clone C-T35 (CABT-ZB313)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes thrombomodulin (TM), a 105 kDa protein
Target	THBD
Immunogen	Eahy 926 cell line
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-T35
Conjugate	Unconjugated
Applications	ELISA(Det) Recommended Working Concentration: ELISA: 0.05-0.5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB309 - CABT-ZB313 CABT-ZB311 - CABT-ZB313
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS

Preservative	None
Storage	Stable at +2-8°C for 12 months
Ship	Wet ice

GENE INFORMATION

Synonyms	THBD; thrombomodulin; TM; THRM; AHUS6; BDCA3; CD141; THPH12; fetomodulin; CD141 antigen
Entrez Gene ID	7056
UniProt ID	P07204