



Mouse Anti-Human CD25 monoclonal antibody, clone C-G3 [FITC] (CABT-ZB319)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes the Interleukin 2 Receptor α chain (IL-2 Ra), a 55 kDa protein
Target	IL2RA
Immunogen	PHA activated T cells
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-G3
Conjugate	FITC
Applications	ELISA(Det), FC, IHC, IP, FA Recommended Working Concentration: ELISA: 0.05-0.5 μ g/ml FC: Use 10 μ l to label 10^6 cells or 100 μ l of whole blood We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB315 - CABT-ZB317; CABT-ZB315 - CABT-ZB319
Format	Purified, Liquid
Concentration	Lot specific
Size	100 tests (1 mL)
Buffer	PBS

Preservative	0.1% Sodium Azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Ship	Wet ice

GENE INFORMATION

Synonyms	IL2RA; interleukin 2 receptor; alpha; p55; CD25; IL2R; TCGFR; IDDM10; interleukin-2 receptor subunit alpha; TAC antigen; IL-2R subunit alpha; IL-2 receptor subunit alpha; HuMax-TAC
Entrez Gene ID	3559
UniProt ID	P01589