



Mouse Anti-Human CD253 monoclonal antibody, clone C-T34 [Biotin] (CABT-ZB323)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes the TRAIL, APO2L antigen, a 32 kDa protein
Target	TNFSF10
Immunogen	Recombinant human TRAIL
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-T34
Conjugate	Biotin
Applications	ELISA(Cap), FC, FA Recommended Working Concentration: ELISA: 0.01-0.1 µg/ml FC: Use 10µl to label 10⁶ cells or 100 µl of whole blood. We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB321 - CABT-ZB326 CABT-ZB323 - CABT-ZB326 CABT-ZB324 - CABT-ZB326
Reconstitution	Reconstitute with 1 mL deionised water.
Format	Purified, Lyophilised
Size	100 tests (1 mL)

Buffer	PBS
Preservative	0.1% Sodium Azide
Storage	Stable at +2-8°C for six months after reconstitution. For longer storage freeze aliquots.
Ship	Wet ice

GENE INFORMATION

Synonyms	TNFSF10; tumor necrosis factor (ligand) superfamily, member 10; TL2; APO2L; CD253; TRAIL; Apo-2L; tumor necrosis factor ligand superfamily member 10; Apo-2 ligand; TNF-related apoptosis inducing ligand TRAIL
Entrez Gene ID	8743
UniProt ID	P50591