



Mouse Anti-Human CD62E monoclonal antibody, clone C-T4 (CABT-ZB337)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes the Endothelial Leucocyte Adhesion Molecule-1 (ELAM-1), E-selectin, a 97 kDa protein
Target	SELE
Immunogen	Activated human endothelial cells
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-T4
Conjugate	Unconjugated
Applications	ELISA(Det) Recommended Working Concentration: ELISA: 0.05-0.5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB335 - CABT-ZB337
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS

Preservative	None
Storage	Stable at +2-8°C for 12 months
Ship	Wet ice

GENE INFORMATION

Synonyms	SELE; selectin E; ELAM; ESEL; CD62E; ELAM1; LECAM2; E-selectin; ELAM-1; endothelial adhesion molecule 1; CD62 antigen-like family member E; endothelial leukocyte adhesion molecule 1; leukocyte endothelial cell adhesion molecule 2; leukocyte-endothelial cell adhesion molecule 2
Entrez Gene ID	6401
UniProt ID	P16581