



Mouse Anti-Human Galectin-9 monoclonal antibody, clone C-L49 (CABT-ZB342)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human Galectin-9
Target	LGALS9
Immunogen	Recombinant human Galectin-9
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-L49
Conjugate	Unconjugated
Applications	ELISA(Det) Recommended Working Concentration: ELISA: 0.05-0.5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB339 - CABT-ZB342
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS
Preservative	None

Storage	Stable at +2-8°C for 12 months
----------------	--------------------------------

Ship	Wet ice
-------------	---------

GENE INFORMATION

Synonyms	LGALS9; lectin, galactoside-binding, soluble, 9; HUAT; LGALS9A; galectin-9; gal-9; ecalectin; galectin 9; tumor antigen HOM-HD-21; urate transporter/channel protein
-----------------	--

Entrez Gene ID	3965
-----------------------	----------------------

UniProt ID	O00182
-------------------	------------------------