



Mouse Anti-Human IFN- γ monoclonal antibody, clone C-H2 (CABT-ZB348)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human IFN-g
Target	IFNG
Immunogen	Recombinant human IFN-g
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-H2
Conjugate	Unconjugated
Applications	ELISA(Det), ELISPOT Recommended Working Concentration: ELISA / ELISpot: 0.05-0.5 μg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB344 - CABT-ZB348 CABT-ZB346 - CABT-ZB348 CABT-ZB347 - CABT-ZB348
Format	Purified, Liquid
Concentration	Lot specific
Size	100 μ g
Buffer	PBS

Preservative	None
Storage	Stable at +2-8°C for 12 months
Ship	Wet ice

GENE INFORMATION

Synonyms	IFNG; interferon, gamma; IFG; IFI; interferon gamma; IFN-gamma; immune interferon
Entrez Gene ID	3458
UniProt ID	P01579