



Mouse Anti-Human IL-10 monoclonal antibody, clone C-U21 (CABT-ZB352)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human IL-10
Target	IL10
Immunogen	Recombinant human IL-10
Isotype	IgG2b, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-U21
Conjugate	Unconjugated
Applications	ELISA(Det), ELISPOT Recommended Working Concentration: ELISA / ELISpot: 0.05-0.5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB350 - CABT-ZB352
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS
Preservative	None

Storage	Stable at +2-8°C for 12 months
----------------	--------------------------------

Ship	Wet ice
-------------	---------

GENE INFORMATION

Synonyms	IL10; interleukin 10; IL10X; interleukin-10; CSIF; IL-10; cytokine synthesis inhibitory factor
-----------------	--

Entrez Gene ID	3586
-----------------------	----------------------

UniProt ID	P22301
-------------------	------------------------