



Mouse Anti-Human IL-12 p35 (IL-12 p70) monoclonal antibody, clone C-U32 (CABT- ZB354)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human IL-12 p35 + p70
Target	IL12
Immunogen	Recombinant human IL-12
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-U32
Conjugate	Unconjugated
Applications	ELISA(Cap), ELISPOT, FA Recommended Working Concentration: ELISA: 0.3-1 µg/ml ELISpot: 5-10 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB354 - CABT-ZB359; CABT-ZB354 - CABT-ZB361
Format	Purified, Liquid
Concentration	Lot specific
Size	1.0 mg

Buffer	PBS
Preservative	None
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Ship	Wet ice

GENE INFORMATION

Synonyms	CLMF; CLMF p35; CLMF p40; CLMF1; CLMF2; Cytotoxic lymphocyte maturation factor 1, p35; Cytotoxic lymphocyte maturation factor 35 kDa subunit; Cytotoxic lymphocyte maturation factor 40 kDa subunit; IL 12 subunit p35; IL 12 subunit p40
Entrez Gene ID	3592
UniProt ID	O60595