



Mouse Anti-Human IL-12 monoclonal antibody, clone C-Q35 [FITC] (CABT-ZB361)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human IL-12 p40 + p70
Target	IL12
Immunogen	Recombinant human IL-12
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-Q35
Conjugate	FITC
Applications	ELISA(Det), ELISPOT Recommended Working Concentration: ELISA / ELISpot: 0.05-0.5 µg/ml FC: It is recommended to establish the optimal dilution according to the system being used. Use 10 µl of antibody to label 5.10⁵ cells We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB354 - CABT-ZB359; CABT-ZB354 - CABT-ZB361
Format	Purified, Liquid
Concentration	Lot specific
Size	100 tests (1 mL)
Buffer	PBS

Preservative	0.1% Sodium Azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Ship	Wet ice

GENE INFORMATION

Synonyms	CLMF; CLMF p35; CLMF p40; CLMF1; CLMF2; Cytotoxic lymphocyte maturation factor 1, p35; Cytotoxic lymphocyte maturation factor 35 kDa subunit; Cytotoxic lymphocyte maturation factor 40 kDa subunit; IL 12 subunit p35; IL 12 subunit p40
Entrez Gene ID	3592
UniProt ID	O60595