



# Mouse Anti-Human IL-13 monoclonal antibody, clone C-C24 (CABT-ZB362)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |                                                                                                                                                                                                                     |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Specificity</b>        | Recognizes both natural and recombinant human IL-13                                                                                                                                                                 |
| <b>Target</b>             | IL13                                                                                                                                                                                                                |
| <b>Immunogen</b>          | Recombinant human IL-13                                                                                                                                                                                             |
| <b>Isotype</b>            | IgG1, κ                                                                                                                                                                                                             |
| <b>Source/Host</b>        | Mouse                                                                                                                                                                                                               |
| <b>Species Reactivity</b> | Human                                                                                                                                                                                                               |
| <b>Clone</b>              | C-C24                                                                                                                                                                                                               |
| <b>Conjugate</b>          | Unconjugated                                                                                                                                                                                                        |
| <b>Applications</b>       | ELISA(Cap), ELISPOT, FA, FC<br>Recommended Working Concentration:<br>ELISA: 0.5-2.5 µg/ml<br>ELISpot: 5-10 µg/ml<br>We recommend the following for sandwich ELISA (Capture - Detection):<br>CABT-ZB362 - CABT-ZB364 |
| <b>Format</b>             | Purified, Liquid                                                                                                                                                                                                    |
| <b>Concentration</b>      | Lot specific                                                                                                                                                                                                        |
| <b>Size</b>               | 1.0 mg                                                                                                                                                                                                              |
| <b>Buffer</b>             | PBS                                                                                                                                                                                                                 |

|                     |                                                                     |
|---------------------|---------------------------------------------------------------------|
| <b>Preservative</b> | None                                                                |
| <b>Storage</b>      | Stable at +2-8°C for 12 months. For longer storage freeze aliquots. |
| <b>Ship</b>         | Wet ice                                                             |

## GENE INFORMATION

|                       |                                                                             |
|-----------------------|-----------------------------------------------------------------------------|
| <b>Synonyms</b>       | IL13; interleukin 13; Il-13; interleukin-13; T-cell activation protein P600 |
| <b>Entrez Gene ID</b> | <a href="#">3596</a>                                                        |
| <b>UniProt ID</b>     | <a href="#">O43644</a>                                                      |