



Mouse Anti-Human IL-13 monoclonal antibody, clone C-Q7 (CABT-ZB364)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human IL-13
Target	IL13
Immunogen	Recombinant human IL-13
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-Q7
Conjugate	Unconjugated
Applications	ELISA(Det), ELISPOT Recommended Working Concentration: ELISA / ELISpot: 0.05-0.5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB362 - CABT-ZB364
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS
Preservative	None

Storage	Stable at +2-8°C for 12 months
Ship	Wet ice

GENE INFORMATION

Synonyms	IL13; interleukin 13; IL-13; interleukin-13; T-cell activation protein P600
Entrez Gene ID	3596
UniProt ID	O43644