



Mouse Anti-Human IL-17F monoclonal antibody, clone C-H57 (CABT-ZB370)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both recombinant and natural human IL-17F
Target	IL17F
Immunogen	Recombinant human IL-17F
Isotype	IgG2b, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-H57
Conjugate	Unconjugated
Applications	ELISA(Cap), ELISPOT Recommended Working Concentration: ELISA: 1-5 µg/ml ELISpot: 5-10 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB370 - CABT-ZB372
Format	Purified, Liquid
Concentration	Lot specific
Size	1.0 mg
Buffer	PBS

Preservative	None
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Ship	Wet ice

GENE INFORMATION

Synonyms	IL17F; interleukin 17F; ML1; ML-1; CANDF6; IL-17F; interleukin-17F; cytokine ML-1
Entrez Gene ID	112744
UniProt ID	Q96PD4