



Mouse Anti-Human IL-17F monoclonal antibody, clone C-G71 (CABT-ZB372)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both recombinant and natural human IL-17F
Target	IL17F
Immunogen	Recombinant human IL-17F
Isotype	IgG2b, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-G71
Conjugate	Unconjugated
Applications	ELISA(Det), ELISPOT We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB370 - CABT-ZB372
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS
Preservative	None
Storage	Stable at +2-8°C for 12 months

Ship

Wet ice

GENE INFORMATION

Synonyms	IL17F; interleukin 17F; ML1; ML-1; CANDF6; IL-17F; interleukin-17F; cytokine ML-1
-----------------	---

Entrez Gene ID	112744
-----------------------	------------------------

UniProt ID	Q96PD4
-------------------	------------------------
