



Mouse Anti-Human IL-23 monoclonal antibody, clone C-A34 (CABT-ZB374)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant IL-23, p19 subunit
Target	IL23A
Immunogen	Recombinant IL-23 (subunit p19)
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-A34
Conjugate	Unconjugated
Applications	ELISA(Cap), IHC, FA Recommended Working Concentration: ELISA: 2.5-5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB374 - CABT-ZB376
Format	Purified, Liquid
Concentration	Lot specific
Size	1.0 mg
Buffer	PBS
Preservative	None

Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
----------------	---

Ship	Wet ice
-------------	---------

GENE INFORMATION

Synonyms	IL23A; interleukin 23; alpha subunit p19; P19; SGRF; IL-23; IL-23A; IL23P19; interleukin-23 subunit alpha; IL-23-A; IL-23p19; IL-23 subunit alpha; interleukin 23 p19 subunit; interleukin-23 subunit p19; JKA3 induced upon T-cell activation; interleukin-six; G-CSF related factor; BI 655066
-----------------	--

Entrez Gene ID	51561
-----------------------	-----------------------

UniProt ID	Q9NPF7
-------------------	------------------------