



Mouse Anti-Human IL-27 monoclonal antibody, clone C-H55 (CABT-ZB378)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human IL-27
Target	IL27
Immunogen	Recombinant human IL-27
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-H55
Conjugate	Unconjugated
Applications	ELISA(Cap), FC Recommended Working Concentration: ELISA: 0.5-2.5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB378 - CABT-ZB380; CABT-ZB378 - CABT-ZB382; CABT-ZB378 - CABT-ZB383
Format	Purified, Liquid
Concentration	Lot specific
Size	1.0 mg
Buffer	PBS
Preservative	None

Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
----------------	---

Ship	Wet ice
-------------	---------

GENE INFORMATION

Synonyms	IL27; interleukin 27; p28; IL30; IL-27; IL27A; IL-27A; IL27p28; interleukin-27 subunit alpha; IL27-A; IL-27-A; interleukin 30; interleukin-30; IL-27 p28 subunit; IL-27 subunit alpha; anti-IL-27
-----------------	---

Entrez Gene ID	246778
-----------------------	------------------------

UniProt ID	Q8NEV9
-------------------	------------------------