



Mouse Anti-Human IL-27 monoclonal antibody, clone C-H50 [PE] (CABT-ZB383)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human IL-27
Target	IL27
Immunogen	Recombinant human IL-27
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-H50
Conjugate	PE
Applications	ELISA(Det) Recommended Working Concentration: ELISA: 0.05-0.5 µg/ml FC: Use 10 µl to label 10⁶ cells or 100 µl of whole blood We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB378 - CABT-ZB380; CABT-ZB378 - CABT-ZB382; CABT-ZB378 - CABT-ZB383
Reconstitution	Reconstitute with 1 mL deionised water.
Format	Purified, Lyophilisated
Size	100 tests (1 mL)
Buffer	PBS

Preservative	0.1% Sodium Azide
Storage	Stable at +2-8°C for six months after reconstitution. DO NOT FREEZE.
Ship	Wet ice

GENE INFORMATION

Synonyms	IL27; interleukin 27; p28; IL30; IL-27; IL27A; IL-27A; IL27p28; interleukin-27 subunit alpha; IL27-A; IL-27-A; interleukin 30; interleukin-30; IL-27 p28 subunit; IL-27 subunit alpha; anti-IL-27
Entrez Gene ID	246778
UniProt ID	Q8NEV9