



Mouse Anti-Human IL-34 monoclonal antibody, clone C-A45 (CABT-ZB385)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes both natural and recombinant human IL-34
Target	IL34
Immunogen	Recombinant human IL-34
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-A45
Conjugate	Unconjugated
Applications	ELISA(Cap) Recommended Working Concentration: ELISA: 1-10 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-ZB385 - CABT-ZB387
Format	Purified, Liquid
Concentration	Lot specific
Size	1.0 mg
Buffer	PBS
Preservative	None

Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
----------------	---

Ship	Wet ice
-------------	---------

GENE INFORMATION

Synonyms	IL34; interleukin 34; IL-34; C16orf77; interleukin-34; anti-IL-34
-----------------	---

Entrez Gene ID	146433
-----------------------	------------------------

UniProt ID	Q6ZMJ4
-------------------	------------------------